
**Flex massagekudde/ Flex massage cushion/
Flex-Massagekissen/ Flex massasjepute/
Flex-massagepuden /FLEX hierontatyyny/
Flexi massaažipatja**

**Manual
User handbook
Handbuch
Brukermanual
Manual
Käsikirja
Juhend**

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info@sableline.se

www.sableline.s

CONGRATULATIONS on your purchase of Flex Massage Cushion from Sableline Healthcare. To get the most out of this product, please read this manual carefully before use.

SAFETY INSTRUCTIONS

Read the following instructions carefully before using FLEX massage cushion, in order to ensure its proper functioning and optimal effect.

This massage cushion complies with the recognised principles of technology and current safety regulations.

Please retain these instructions for use.

- All repairs must be carried out by qualified persons. Faulty use and unauthorised repairs of the product are for safety reasons prohibited and will render the warranty invalid.
- Never touch the electrical plug with wet hands.
- Do not allow FLEX massage cushion to come into contact with water or direct sunlight.
- Make a point of disconnecting the massage cushion from the mains power socket when it is going to be out of use for an extended period.
- Do not use damaged plugs or cord. A damaged plug must be replaced by the manufacturer or service representative (warranty period 2 years), after which it can be purchased from the manufacturer (see warranty). In the event of a malfunction, disconnect the FLEX massage cushion electrical adapter immediately. We are unable to assume liability for damage resulting from incorrect handling.
- In order to avoid strained muscles and nerves, the recommended period of use should not be exceeded.
- In order to prevent unnecessary damage to FLEX massage cushion, never use sharp or pointed objects in conjunction with FLEX massage cushion.
- All massage with FLEX massage cushion must be avoided during pregnancy and/or if pain is experienced in the area to be massaged which is associated with: a disease of the head, a recent injury, thrombotic illnesses, inflammations and swellings of any kind, as well as cancer. With other serious diseases, consultation with your doctor is recommended.
- Do not use FLEX massage cushion while you are driving.
- If you depend on electronic aids such as a pacemaker, consult your doctor before use.

INSIDE THE PACKAGE

- MC – G5 FLEX massage cushion
- Mains power adapter
- 12 V car adapter
- Active straps

TECHNICAL DATA

- MC-G5 FLEX massage cushion.
- Article No. 14-600-00.
- Timer 15 min.
- Input 100-240 V 50/60 Hz 0.8 A
- Output 12 V 2500mA
- Rated power 30W
- Certificate

INSTRUCTIONS FOR USE

MC-G5 FLEX MESSAGE CUSHION CAN BE USED TO BENEFICIAL EFFECT ON ALL AREAS OF THE BODY

Neck – Shoulders – Upper back – Lower back – Thighs – Calves – Feet – Hands – Arms

FLEX, your practical and flexible massage cushion, will gently loosen up your muscles and quickly increase your general feeling of well-being without the need for you to get up from your sofa or armchair. Create an oasis of wellness in the comfort of your own home where all the stresses of mind and body are resolved. The four ergonomically designed massage heads in FLEX massage cushion provide a massage of calming circular movements on the desired area of the body. With the aid of the heating function, the muscles are gently warmed up and within several minutes a feeling of well-being is generated which permeates the whole body. Regular massage with FLEX massage cushion will also stimulate the nerves and increase the metabolism. It stimulates the circulation, releases stubborn tensions and has a positive impact on blockages or disorders of internal organs. Shiatsu massage always has a positive effect, regardless of whether it is intended to prevent problems or simply for relaxation and activation of life energy. Regular and well-applied massage improves the circulation and provides alleviation of pain in the neck and shoulders, and further strengthens the large muscle groups in the back. Massage also stimulates the healthy functioning of the nerves and increases the body's metabolism. Massage of reflex points on the bottom of the feet generates a positive feeling of harmony throughout the body.

SHIATSU MASSAGE releases deep-seated tensions, strengthens the back muscles and activates lymph drainage and the circulation.

HEAT enhances the penetrative effect and softens up the musculature.

VIBRATION is effective for chronic pain and increases the circulation.

Four massage heads perform circular movements to provide a very comfortable massaging effect on the neck, shoulders, upper back, lower back, thighs, calves, hands and feet. FLEX massage cushion massages in two directions, clockwise and anti-clockwise, and further incorporates infrared heating and vibration. FLEX massage cushion can be connected to a 220 V or 12 V (adapter for 12 V included) power source, and is therefore also usable in cars, boats, caravans and motor caravans.

When massaging the neck, upper arms, upper back and lower back it is best to use the active straps included in the package. The active straps ensure optimal positioning of FLEX massage cushion on the neck, shoulders, upper arms, upper back and lower back, and facilitate adjustment of FLEX, making it easy for you to reach those areas that need massaging. Another advantage of the active straps is that they make it easy to regulate the pressure and to move FLEX massage cushion sideways in order to manoeuvre it into the exact position you want. This means it is easier for you to massage your neck, upper arms and shoulders without requiring the support of a high back on your sofa or armchair (see illustrations).

Underarms

Hands

Thighs (front)

Thighs (sides)

Thighs (back)

Feet

HOW TO MASSAGE

Underarms: lay FLEX massage cushion on the armrest of the sofa or armchair and then massage one arm at a time.

Hands: lay the Flex massage cushion on the knee and lay your hands on FLEX.

Thighs (back): easiest to massage by laying FLEX massage cushion on the seat of the sofa or armchair and then massaging one thigh at a time. You can of course lie down on your bed or sofa to do the same thing.

Thighs (front): lay FLEX massage cushion on top of your thigh and press down on it to achieve the desired effect.

Thighs (sides): lay FLEX massage cushion on your sofa or armchair so that it is supported by the arm rest.

Hip: lay FLEX massage cushion diagonally in the corner of your sofa or armchair and press your hip against it.

Calves: these are easiest to massage by using a footstool, or again by lying on your bed or sofa.

Feet: lay FLEX massage cushion in front of your sofa or armchair and massage one or both feet at the same time (see illustrations).

ACTIVE STRAPS

The active straps are easily connected to FLEX massage cushion with the silver-coloured loops which are located on the sides of FLEX massage cushion. The length of the active straps is easily adjusted with the aid of the Velcro® fastening on the inside of the straps (see illustrations).

- Do not use FLEX massage cushion on the body for longer than 30 minutes.
- After massage is completed, wait for approximately 2 hours before re-using.
- Initially, we recommend that you do not use FLEX massage cushion for more than 3–5 minutes on each part of the body, as doing so could result in training pain.
- Avoid massaging the spine directly.
- Medical recommendations: unless your doctor or orthopaedist has expressly told you otherwise, then it is safe for you to use FLEX massage cushion. Massage must never feel unpleasant. Excessive massage could be harmful to your body. Pregnant women, persons who are ill and small children should not use the cushion. For persons with a spinal injury or who have recently undergone an operation, massage could be directly harmful and should therefore be avoided. If you are uncertain, always consult your doctor.

CONTROL PANEL

Pressing the **ON/OFF** button will start FLEX massage cushion. Once started, the infrared heating function will automatically be activated. If you want to stop the heating function, press the ON/OFF button once more. A timer will shut down FLEX massage cushion after 15 minutes – if further massage is required, press the ON/OFF button again. If you want to stop massaging before the 15 minutes are up, hold the ON/OFF button down for approximately 5-10 seconds.

Pressing **VIBRATION** activates the vibration function in FLEX massage cushion – this can only be used in conjunction with Shiatsu massage.

TROUBLESHOOTING

Problem

FLEX massage cushion does not function even though the mains power adapter is connected.

Justering

Check that the mains power adapter is connected.

Check whether FLEX massage cushion functions with the 12 V adapter, in order to establish whether the fault lies with the mains power adapter.

FLEX massage cushion can become overheated when it is turned off automatically. Allow to cool for approximately 30-40 minutes and then it will resume normal functioning.

SAFETY INSTRUCTIONS

MC – G5 FLEX massage cushion has an approved CE marking and complies with the recognised principles of technology and current safety regulations. Bearing in mind the relevant safety aspects, you should observe the following instructions:

- Separate FLEX massage cushion from its adapter after use
- Do not pull on the cord when removing the mains power adapter from the wall socket.
- Use only approved electrical sockets.
- When caring for your FLEX massage cushion use just a soft towel or a mild detergent.
- Using the zip fastener on the side of FLEX massage cushion you can open up the cushion and clean inside.
- Never use solvents.

El-Kretsen →

WARRANTY CONDITIONS

WARRANTY

Provided that MC-G5 FLEX massage cushion is not subjected to direct damage and that the care instructions are followed, a 2-year warranty is provided on both motor and material.

Date of purchase: _____

For easier and faster handling of warranty claims, please use the e-mail address: info@sableline.se.

Describe the fault and give your name, address and telephone number, and attach a copy of the warranty certificate or purchase receipt.

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Sableline Sweden AB